## INSTRUCTION FOR USE

- en Instruction for use
- bg Инструкции за употреба
- cs Návod k použití
- da Brugsanvisning
- de Gebrauchsanweisung
- es Instrucciones de uso
- fi Käyttöohjeet
- fr Mode d'emploi
- el Οδηγίες χρήσης
- hu Használati utasítás
- it Istruzioni per l'uso
- lt Naudojimo instrukcija
- nl Gebruiksaanwijzing
- no Bruksanvisning
- pl Instrukcja użytkowania
- pt Instruções de utilização
- ro Instrucțiuni de utilizare
- sk Návod na použitie
- sv Bruksanvisning
- tr Kullanım Kılavuzu
- uk Інструкція з експлуатації

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Cold temperature indicator: Do not use product if indicator bulb is purple

· Индикатор за прекалено ниска температура: Продуктът да не се използва, ако индикаторната крушка е в лилаво · Indikátor nízkej teploty: Produkt nepoužívejte, pokud je žárovka indikátoru nachová · Indikator for lav temperatur: Produktet må ikke anvendes, hvis indikatorpæren er lilla · Kälteindikator: Das Produkt nicht verwenden, wenn die IndikatorKugel violett ist · Indicador de temperatura baja: No utilizar el producto si el indicador es de color violeta · Kylmän lämpötilan osoitin: tuotetta ei saa käyttää, jos osoitin on purppuranvärinen · Témoin de température basse : ne pas utiliser le produit si le témoin est violet · Ένδειξη χαμηλής θερμοκρασίας: Μην χρησιμοποιείτε το προϊόν εάν ο ενδεικτικός λαμπτήρας είναι μοβ · Hideghőmérséklet-jelző: Ne használja a terméket, ha a jelző lila. · Indicatore di temperatura fredda: non utilizzare il prodotto se l'indicatore è viola · Šaltos temperatūros indikatorius: Nenaudokite gaminio, jei indikatorius lemputė violetinė · Koudetemperatuurindicator: Het product niet gebruiken als het bolletje van de indicator paars is · Indikator for kald temperatur: Bruk ikke produktet hvis indikatorens blære er lilla · Wskaźnik niskiej temperatury: nie używać produktu, jeśli żarówka wskaźnika jest fioletowa · Indicador de temperatura fría: não utilizar o produto se a lâmpada do indicador estiver roxa · Indicator pentru temperatură joasă: nu utilizați produsul dacă indicatorul este violet · Indikátor nízkej teploty: Nepoužívajte produkt, ak je indikátor fialový · Indikator för kall temperatur: Använd inte produkten om indikatorn är lila · Düşük sıcaklık göstergesi: Göstergenin lambası morsa ürünü kullanmayın · Индикатор холодных температур: Не используйте вириб, якщо колба індикатора фіолетова

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
**Please read all instructions carefully before using the product. Please comply with all warnings and precautions included in these instructions. Non-compliance may result in serious complications.**

**Implantation of the ALLEGRA Transcatheter Heart Valve may only be performed by physicians who have received appropriate training.**

#### 1. DEVICE DESCRIPTION

The ALLEGRA Transcatheter Aortic Valve Implantation System transfemoral approach (ALLEGRA TAVI System TF) consists of the ALLEGRA Transcatheter Heart Valve (REF see chapter 1.1) and the ALLEGRA Delivery System TF (REF DSL- AO18G1RE1150) including the ALLEGRA Loading System TF.

##### 1.1. ALLEGRA TRANSCATHETER HEART VALVE (BIOPROSTHESIS)

The ALLEGRA Transcatheter Heart Valve (THV) is intended to be implanted within degenerative calcified aortic valves or in failing surgical aortic valve bioprostheses in a minimal invasive transcatheter implantation technique.

The supra-annular valve has a tri-leaflet design and is constructed with six individual bovine pericardium parts, three parts for the skirt and three leaflets sewn to the skirt in a semi-lunar fashion to form the valve occluder. The pericardial inner skirt creates a 12 mm sealing against the native aortic annulus and thereby mitigating the risk of paravalvular leakage. The bovine pericardium used for the THV is treated to reduce calcification potential. The stent is a nitinol laser-cut stent, with good radiopaque visibility and 6 additional gold radiopaque markers, where the valve is sutured in. The table below indicates the available sizes of the THV with catalog numbers and the use in regard of the native aortic annulus diameter of the patient and the corresponding true inner diameter of a degenerated surgical aortic valve bioprosthesis.

<b>ALLEGRA Transcatheter Heart Valve (THV)</b>	<b>ALLEGRA 23</b>	<b>ALLEGRA 27</b>	<b>ALLEGRA 31</b>
<b>Catalog / Reference number</b>	THV-AO0023G1RE	THV-AO0027G1RE	THV-AO0031G1RE
<b>Aortic annulus diameter [mm]</b>	19 - 22	22 - 25	25 - 28
<b>Inner diameter aortic valve bioprosthesis [mm]</b>	16.5 - 21.5	22 - 24.5	25 - 28

##### 1.2. ALLEGRA DELIVERY SYSTEM TF (DELIVERY SYSTEM)

The ALLEGRA Delivery System TF is intended for transfemoral retrograde implantation of the ALLEGRA Transcatheter Heart Valve (THV).

The delivery system is compatible with a 0.035 inch (0.889 mm) guidewire and has to be inserted to the vessel through a separate introducer sheath. The catheter shaft has a distal short cartridge of 18 Fr, which accommodates the THV. Behind the cartridge section, the catheter shaft reduces to 15 Fr, which gives the necessary flexibility, stability and pushability. The usable length of the catheter shaft (1150 mm) allows the implantation site to be reached in the target patient population. The atraumatic catheter tip and the loading connector are radiopaque, and the inflow cover tip and the cartridge of the catheter shaft have radiopaque rings. The delivery system incorporates a proprietary technology Permaflow which allows an occlusion free deployment and

permits repositioning and retrieval of the THV in Permaflow. In Permaflow the valve is unsheathed to 80% leaving the inflow and outflow retained within the delivery catheter. In this phase, the valve can be fully recaptured inside the cartridge of the catheter by advancing the retaining catheter over the exposed valve. For opening the inflow of the THV the catheter tip is advanced until the inflow portion of the valve is fully released. At this stage the valve is fully functional and >90% of the stent frame is released from the catheter. If the valve position is not acceptable, the valve can still be recaptured and retrieved from the patient using the procedure described in chapter 7.13.

Model	Catalog / Reference number	Compatible with ALLEGRA THV
ALLEGRA Delivery System TF	DSL- AO18G1RE1150	ALLEGRA 23, ALLEGRA 27, ALLEGRA 31

### 1.3. ALLEGRA LOADING SYSTEM TF (LOADING SYSTEM)

The ALLEGRA Loading System TF is intended to load the ALLEGRA Transcatheter Heart Valve into the cartridge of the ALLEGRA Delivery System TF.

The loading system consists of a loading tool, a counter funnel, a loading pin to orientate the pericardium inside the stent, and a loading clamp to hold the tip release mechanism during loading. The loading system is included in the packaging of the delivery system.

## 2. INDICATION

The ALLEGRA Transcatheter Aortic Valve Implantation System TF is indicated for the treatment of severe calcified aortic valve stenosis in high risk patients with elevated surgical risk or in patients with a symptomatic degeneration of an aortic valve bioprosthesis.

### 2.1. TARGET PATIENT GROUP

The device was studied in high risk patients (STS-PROM / EuroSCORE II >8%) with symptomatic, severe native aortic stenosis or in symptomatic patients with failing surgical aortic bioprosthetic valves (valve-in-valve procedure) who were considered at high risk for re-do surgery, as assessed by a multidisciplinary heart team. The degeneration of the surgical bioprosthetic valve in the valve-in-valve patient population resulted from valve stenosis, insufficiency or the combination of both. The target patient group shall present:

- Symptomatic degeneration of severe calcified aortic heart valve or failing surgical aortic bioprosthetic valve showing echocardiographically mean aortic gradient  $\geq 40$ mmHg or peak jet velocity  $\geq 4.0$ m/s or AVA  $\leq 1.0$  cm<sup>2</sup> (or  $\leq 0.6$  cm<sup>2</sup>/m<sup>2</sup>) or symptomatic patients with severe bioprosthetic valve insufficiency.
- Current heart valve replacement guidelines (e.g. Guideline of European Society of Cardiology) recommend bioprosthetic valves for patients  $\geq 75$  years of age or above.
- Anatomic requirements according to device specifications in chapter 1.1. and 1.2.

## 3. CONTRAINDICATION

The ALLEGRA Transcatheter Aortic Valve Implantation System TF must not be used if the implanting physician believes its implantation would be contrary to the best interest of the patient.

### 3.1. ABSOLUTE CONTRAINDICATIONS

The ALLEGRA TAVI System TF is contraindicated in the presence of any of the following conditions:

- Unicuspid or bicuspid aortic valve.
- Non-calcified aortic stenosis.

- Combined aortic valve disease with predominant aortic regurgitation > 3.
- Distance between aortic valve basal plane and the orifice of the lowest coronary artery < 8 mm.
- Echocardiographic evidence of intracardiac thrombus or vegetation.
- Significant aortic disease such as severe obstructive calcification or marked tortuosity or kinking which would preclude a safe advancement of the ALLEGRA TAVI System TF.
- Iliofemoral vessel conditions such as severe obstructive calcification, severe tortuosity or kinking that would preclude safe placement of an 18 Fr introducer sheath or make endovascular access to the aortic valve impossible.
- Severe ventricular dysfunction with LVEF less than 20%.
- Evidence of active endocarditis or other acute infections.
- Renal failure requiring continuous renal replacement therapy.
- Known hypersensitivity to contrast media, which cannot be adequately pre-medicated or contraindication to anticoagulant or anti-platelet medication (e.g. aspirin, heparin, or clopidogrel) or to Nitinol (titanium or nickel) alloy or to bovine tissue.

### 3.2. RELATIVE CONTRAINDICATIONS

- Evidence of an acute myocardial infarction within the past 30 days.
- Cerebral vascular accident (CVA) within past 6 months.
- Evidence of active peptic ulcer or upper gastrointestinal bleeding within past 90 days.
- Severe mitral regurgitation.

### 4. WARNINGS

- This procedure should only be performed where emergency aortic valve surgery can be performed promptly.
- Correct sizing of the THV is essential to prevent paravalvular leakage, migration or annular rupture. The THV is intended for use in candidates with a native aortic annulus size ranging from 19 mm to 28 mm or with a failing surgical aortic valve bioprosthesis with a true inner diameter ranging from 16.5 mm to 28 mm.
- Correct implantation height of the THV is essential to prevent interferences with the mitral leaflet, coronary obstruction, and poor performance of the THV.
- To avoid occlusion of the coronary ostia, a distance between aortic valve basal plane and the orifice of the lowest coronary artery < 8 mm must be avoided in order to have a safety margin and not to cover the ostia with the skirt of the THV. The upper limit of the skirt is visible through 6 radiopaque markers.
- Do not resterilize or re-use the THV, the delivery system or the loading system. Resterilization or re-use may compromise the structural integrity of the products and/or create a risk of contamination of the device, which could result in patient injury, illness or death.
- Do not use the products if the use by date has expired
- Accelerated deterioration due to calcific degeneration of the THV (as with any glutaraldehyde cross-linked bioprosthesis) may occur in patients with an altered calcium metabolism.
- Overall durability, especially long-term, has not been established for the THV. Careful and continuous medical follow-up is advised so that THV-related complications can be diagnosed and properly managed.
- It is recommended that all THV recipients be prophylactically treated for endocarditis to minimize the possibility of prosthetic valve infection.

- THV recipients should be maintained on anticoagulant therapy, except where contraindicated, as determined by their physician.
- Do not freeze. Always store the THV in a cool, dry place. Any THV that has been frozen must not be used for human implantation.
- The 18 Fr introducer sheath from Cook Medical must not be used because it may cause damage to the tip of the cartridge.
- Do not perform an implantation without an introducer sheath or a guidewire.
- Physicians without previous experience using the ALLEGRA TAVI System TF shall be trained according to the NVT Physician Training Plan.

## 5. PRECAUTIONS

### 5.1. PRECAUTIONS PRIOR TO USE

- Before, during and post-implant handling of the ALLEGRA TAVI System TF shall be in accordance with instructions provided in chapter 7 “Directions for use”.
- Before using the products, check the temperature indicators, the ETO indicator, the compliance of the series or lot number, and the correct size of the implant by using the four-eye principle.
- Do not use the THV if the warm temperature indicator is black or the cold temperature indicator’s bulb is purple.
- Use the delivery system only if the ETO indicator on the inner sterile bag is green.
- Do not use if the integrity of the sterile packaging has been compromised.
- In the event of any non-matching information, do not use the product.
- Do not use the THV if the product jar is damaged or if the seal is broken or the outside of the product jar is wet.
- Do not use the THV if the level of the liquid inside the product jar is below the 120 ml mark.
- Do not use the THV if there are any signs of damage to the tissue or frame.
- Do not use the delivery system if there are any signs of damage.

### 5.2. PRECAUTIONS DURING USE

- The outside of the product jar is not sterile and must not enter the sterile field.
- Follow aseptic techniques during device preparation and implantation.
- Any THV that is dropped, damaged, or mishandled in any way must not be used for implantation.
- Do not handle or manipulate the THV with sharp or pointed objects.
- Do not touch or squeeze the THV during rinsing. No other objects should be placed in the rinsing bowls.
- Make sure not to rinse the THV in an already used rinsing solution.
- Do not load the THV without protection wire.
- Do not use the delivery system if flushing of the system is not possible.
- Avoid too many crimping cycles during loading.
- Make sure that the THV touches the catheter tip, otherwise the reference with the radiopaque ring will be wrong.
- Check if the safety locker is locked and cannot be moved.
- If you feel any resistance during loading, the loading procedure must be stopped and the cause of the resistance evaluated.
- Do not use any tool other than the loading pin, in order to prevent damage to the pericardium.
- Do not use any tool other than the loading pin, to press the struts inside.
- Make sure that the suture or particles from the heart valve holder do not remain inside the THV.

- Make sure that no parts of the tip funnel remain inside or on the catheter.
- Keep the THV hydrated and de-aired by flushing the port of the catheter shaft until ready for implantation.
- In case of an implantation in a failing surgical bioprosthesis, carefully check the patient's anatomy regarding low coronary artery location, especially in combination with shallow sinuses.
- In case of an implantation in a failing surgical bioprosthesis, ensure that there is no partially detached valve leaflet, which can occlude a coronary ostium in the aortic position.
- In case of an implantation in a failing surgical bioprosthesis, the risk of coronary occlusion is higher for surgical bioprostheses with external leaflets.
- In case of an implantation in a failing surgical bioprosthesis, ensure that the guidewire is in the central lumen of the failing surgical bioprosthesis because paravalvular placement or placement through a defect in the leaflet or the suture ring of the failing surgical bioprosthesis will lead to complications.
- Ensure the diameter of the access vessels is adequate.
- Do not kink the catheter shaft.
- Make sure that if you pre-bend the guidewire it will not be damaged.
- Balloon aortic valvuloplasty (BAV) of the native aortic valve is required prior to delivery system insertion to facilitate placement of the THV.
- The balloon valvuloplasty catheter should be sized at least 1 to 2 mm smaller than the annular diameter. For the ALLEGRA 31 at least a 24 mm balloon is recommended.
- Balloon valvuloplasty should only be performed when the THV is loaded into the delivery system and ready for implantation.
- In case of failing surgical bioprostheses with structural damage, balloon valvuloplasty should be performed with extreme caution.
- Before introducing the delivery system into the patient, check visually the correct loading of the THV.
- If you feel any resistance introducing the delivery system, the procedure must be stopped and the cause of the resistance has to be evaluated.
- Make sure that the catheter at the integrated sheath is stable, once positioned in correct annular position in order to avoid any change.
- Re-sheath the THV only in bailout situations. Do not re-implant the THV after re-sheathing.
- Make sure that the guidewire does not hinder the withdrawal of the catheter tip.
- Before removal of the delivery system, check if the T-bars are free from the catheter shaft after release of the THV.
- The catheter closure is recommended to be done in a straight portion of the descending aorta.
- During catheter closure do not force the catheter against the tip to avoid overlapping. If this occurs the diameter of the tip will be increased by the overlapping catheter and makes it difficult to cross the introducer sheath.
- Proceed with caution if it is necessary to post-dilate the THV in the case of unacceptable paravalvular leakage. Post-dilatation could damage the device integrity or cause migration of the THV.
- Glutaraldehyde may cause irritation of the skin, eyes, nose and throat. Avoid prolonged or repeated exposure or inhalation of the solution. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with water. In the event of contact with eyes, seek immediate medical attention.

## 6. POTENTIAL COMPLICATIONS, RISKS AND SIDE EFFECTS

Potential complications, risks and side effects that may be associated with the use of the ALLEGRA TAVI System TF include, but may not be limited to:

- Acute coronary occlusion
- Acute kidney injury
- Acute myocardial infarction
- Acute renal failure
- Allergic reactions/intolerances (e.g. to contrast media)
- Aortic root injury (e.g. dissection, perforation)
- Arrhythmias including ventricular tachycardia or fibrillation extending to cardiac arrest
- Atrioventricular conduction disorders (e.g. AV-block, LBBB) which may require a permanent pacemaker
- Bleeding (hemorrhage)
- Cardiac tamponade
- Cardiogenic shock
- Cerebrovascular events such as TIA, Stroke
- Death
- Delivery system malfunction requiring the exchange of the device
- Device embolization / migration
- Emergent cardiac surgery (e.g. coronary artery bypass or heart valve replacement) or additional transcatheter intervention (e.g. second valve implantation or change to another device)
- Endocarditis
- Exacerbation of heart failure
- Hemolysis
- Hemorrhage requiring transfusion
- Hypertension or hypotension
- Infection
- Mitral valve injury
- Non-structural prosthetic valve dysfunction: paravalvular or/and central regurgitation, stenosis caused by mis-sizing and/or patient prosthesis mismatch
- Prosthetic valve thrombosis
- Sepsis
- Structural prosthetic valve damage (e.g. cusp tear, suture line disruption, stent fracture, calcification)
- Thromboembolism
- Vascular injury (e.g. dissection, perforation)

## 7. DIRECTIONS FOR USE

### 7.1. PHYSICIAN TRAINING

The implanting physician should be experienced in port access, balloon aortic valvuloplasty, catheterization techniques, the THV implantation procedure, and trained on the use of the THV, the delivery system and the loading system.

### 7.2. REQUIRED EQUIPMENT

Standard cardiac catheterization lab equipment is required.

### 7.3. REQUIRED MATERIAL

The materials listed below are not delivered with the ALLEGRA TAVI System TF:

- 3 rinsing bowls and 1 loading bowl (capacity approx. 1 l / bowl)
- 1 small bowl for heparinized saline flushing solution (capacity approx. 100 ml)
- 2 liter sterile 0.9% saline solution
- Approx. 50 ml heparinized sterile 0.9% saline solution
- Sterile ice cubes
- Guidewire, 0.035 inch (0.889 mm) (recommended: extra-stiff guidewire of 300 cm length)
- Percutaneous transluminal valvuloplasty catheter in various sizes
- Compatible 18 Fr introducer sheath
- Femoral closure device
- Pigtail catheter
- 10 ml straight tip syringe (for catheter flushing)
- Standard procedure supplies

**CAUTION:** The 18 Fr introducer sheath from Cook MEDICAL must not be used because it may cause damage to the tip of the cartridge.

### 7.4. INSPECTION PRIOR TO USE

**NOTE:** Before using the products, check the temperature indicators, the ETO indicator, the compliance of the series or lot number, and the correct size of the implant by using the four-eye principle.

- Carefully inspect the package before opening.

**CAUTION:** Do not use if the integrity of the sterile package has been compromised. **CAUTION:** Do not use the THV if one of the temperature indicators has been activated.

- Remove the product from the protective package and visually check that it is free of defects.

**CAUTION:** Do not use if any defects are detected.

**CAUTION:** Use the delivery system only if the ETO indicator on the inner sterile bag is green.

### 7.5. THV HANDLING AND PREPARATION

**NOTE:** Follow aseptic techniques during device preparation and implantation.

Procedure	
1	Remove the product jar from the carton box. Verify that all carton box and product jar labels match with respect to THV model, size and serial number. <b>CAUTION:</b> In the event of any non-matching information, do not use the product.
2	Carefully inspect the entire product jar and seal for damage. <b>CAUTION:</b> Do not use the THV if the product jar is damaged or if the seal is broken or the outside of the product jar is wet. <b>CAUTION:</b> Do not use the THV if the level of the liquid inside the product jar is below the 120 ml mark.
3	Before opening the product jar, prepare three sterile bowls, each containing at least 500 ml of sterile 0.9% saline solution at ambient temperature (15°C to 25°C) to thoroughly rinse the glutaraldehyde sterilant from the THV, and one sterile bowl containing at least 500 ml of cold, sterile 0.9% saline solution (0°C to 8°C) for the loading procedure.

Procedure	
4	Break the seal from the product jar and unscrew the screw cap. The content of the product jar is sterile and must be handled aseptically to prevent contamination. <b>CAUTION:</b> The outside of the product jar is not sterile and must not enter the sterile field.
5	Open and hold the product jar, while someone inside the sterile field takes out the transcatheter heart valve holder from the product jar. The THV must be inspected for any signs of damage to the frame or tissue. <b>CAUTION:</b> Do not use the THV if there are any signs of damage to the tissue or frame. <b>CAUTION:</b> Any THV that is dropped, damaged, or mishandled in any way must not be used for implantation.

### 7.6. THV RINSING PROCEDURE

**NOTE:** The whole rinsing and loading procedure requires approx. 10 minutes; 6 minutes for rinsing (3 x 2 minutes) and approx. 4 minutes for loading.

Procedure	
1	The THV needs to be rinsed to be free from residues of glutaraldehyde. Place the THV in the first bowl of sterile 0.9% saline solution. Be sure the solution completely covers the THV and holder. <b>CAUTION:</b> Do not handle or manipulate the THV with sharp or pointed objects.
2	With the THV and holder submerged, gently agitate the THV back and forth by grasping the holder for a minimum of two minutes. <b>CAUTION:</b> Do not touch or squeeze the THV during rinsing. No other objects should be placed in the rinsing bowls.
3	Repeat the two-minute rinse in the second and third bowl for a total rinse time of at least six minutes. <b>NOTE:</b> Make sure not to rinse the THV in an already used rinsing solution.
4	Leave the THV in the final rinsing solution until needed to prevent the tissue from drying. <b>CAUTION:</b> The THV should be kept hydrated throughout the rest of the preparation procedure to prevent the tissue from drying.

### 7.7. DELIVERY SYSTEM HANDLING AND PREPARATION

Procedure	
1	Carefully inspect the entire sterile bags for damage. <b>CAUTION:</b> Do not use the delivery system if the sterile bags are damaged. <b>CAUTION:</b> Use the delivery system only if the ETO indicator on the inner sterile bag is green.
2	Remove the double sterile bag from the carton box. Verify that all carton box and sterile bag labels match with respect to model, size and lot number. <b>CAUTION:</b> In the event of any non-matching information, do not use the product.

Procedure	
3	<p>Open the outer sterile bag and ensure that only the inner sterile bag is introduced into the sterile field. Remove the delivery system from the inner sterile bag. Remove the shrink seal from the handle and inspect the product for any signs of damage.</p> <p><b>CAUTION:</b> Do not use the delivery system if there are any signs of damage.</p>

### 7.8. THV LOADING PROCEDURE

**CAUTION:** Do not load the THV without a protection wire.

**CAUTION:** Do not use the delivery system if flushing of the system is not possible.

**CAUTION:** If you feel any resistance, the loading procedure must be stopped and the cause of the resistance evaluated.

**CAUTION:** Avoid too many crimping cycles during loading.

Procedure	
1	<p>Unlock the safety locker by pressing both safety locker triggers inwards and move the safety locker backwards until it touches the end of the safety wings.</p>
2	<p>Move the switch on the handle from the position P (Park) to N (Neutral) and move the handle backward until the loading connector is exposed by the retracted catheter shaft.</p>
3	<p>Remove the loading plate from the loading guides. Place the THV on the loading tool, with the T-Bars pointing towards inside the loading funnel. Pass the guiding pins of the loading tool through the loading guides of the heart valve holder. Mount the loading plate on the loading guides.</p> <p><b>NOTE:</b> For easier handling, bend the jar fixation of the valve holder 90° to the side.</p> <p><b>NOTE:</b> Check if the Centering hub is held in position by the spacer.</p>
4	<p>Introduce the tip of the catheter into the funnel hole of the loading tool. Press the THV with the loading plate against the funnel to carefully align the three T-Bars to the grooves of the loading connector and guide them into place. The T-Bars must stay in the grooves.</p> <p><b>NOTE:</b> Make sure not to manipulate the leaflets of the THV.</p>
5	<p>Press both safety locker triggers inwards and move the safety locker with the handle gently forward until it stops and reaches the end position. Release the safety locker triggers. The T-Bars of the stent are secured by the cartridge.</p> <p><b>NOTE:</b> Check if the safety locker is locked and cannot be moved.</p>
6	<p>Remove the loading plate and move the remaining loading tool into the direction of the handle. It is no longer needed and must be removed after the loading process is finished. Cut the three sutures of the heart valve holder and remove the holder.</p>
7	<p>Using the loading pin, push the pericardium, which protrudes outside the rhombs, gently back inside the stentframe.</p> <p><b>CAUTION:</b> Do not use any tool other than the loading pin, in order to prevent damage to the pericardium.</p>

	<b>Procedure</b>
8	<p>Drive the THV inside the cartridge by moving the handle until the flushing port of the catheter shaft is accessible (approximately 1 cm). To secure the reached position push the lever of the handle down and move the switch from the position N (Neutral) to the position P (Park).</p> <p><b>CAUTION:</b> Make sure that the suture or particles from the heart valve holder do not remain on the THV.</p>
9	<p>Cool the THV inside a bowl with cold sterile 0.9% saline solution and continue the loading procedure. Take the counter funnel and crimp the THV inside until the tip of the protection wire touches the counter funnel. Then remove the counter funnel.</p>
10	<p>Push the tip actuator forward to bring the tip in front of the THV. Lock the tip actuator to the tip grip using the loading clamp.</p> <p><b>NOTE:</b> Check if the safety wings came out. If they did not, verify the correct position of the safety locker.</p> <p><b>NOTE:</b> The Centering hub must remain inside the THV and inside the cartridge. The Centering hub should not be visible.</p> <p><b>NOTE:</b> The loading spacer must align with the tip funnel.</p>
11	<p>Before moving on with the loading procedure, a pre-flushing of the catheter shaft must be performed while the THV is partially loaded. For this purpose, use a 10 ml syringe filled with heparinized, sterile 0.9% saline solution and rinse the catheter shaft using the flushing port which is between the handle and the safety locker.</p> <p><b>NOTE:</b> To ensure that the catheter shaft is safely flushed, the tip of the catheter must be held in a higher position than the handle. The pre-flushing will facilitate the final flushing process.</p>
12	<p>Use the loading pin to gently reposition any pericardium, which may be protruding through the rhombs of the stent frame.</p> <p><b>CAUTION:</b> Do not use any tool other than the loading pin, in order to prevent damage to the pericardium.</p>
13	<p>Remove the spacer.</p>
14	<p>Take the protection wire out of the delivery system.</p>
15	<p>Place the THV inside a bowl with cold sterile 0.9% saline solution and crimp it completely into the cylindrical part of the counter funnel. Do not remove the counter funnel. Move the switch on the handle from the position P (Park) to N (Neutral) and move the catheter shaft with the handle slowly forward until only two millimeters of the THV remain visible. Remove the counter funnel. Press three evenly distributed struts inside the other struts with the loading pin.</p> <p><b>CAUTION:</b> Do not use any tool other than the loading pin, to press the struts inside.</p>

	<b>Procedure</b>
16	<p>Take off the loading clamp and carefully release the tip actuator until the tip funnel of the catheter tip encloses the THV. Move the catheter shaft slightly backwards with the handle and press the THV into the inflow cover tip until the struts are at the level of the radiopaque ring.</p> <p><b>NOTE:</b> Prevent holding the tip while releasing the tip actuator. <b>NOTE:</b> Make sure that the tip funnel is concentrically on the stent. <b>NOTE:</b> Avoid torquing or tilting the tip, only press the tip and bring it over the THV.</p>
17	<p>Rotate the delivery system until the fixation of the tip funnel is facing upwards. Press on the metal fixation of the tip funnel until the tip funnel breaks underneath. Put no direct pressure onto the tip.</p> <p><b>NOTE:</b> Make sure not to pull the fixation to remove the tip funnel. <b>CAUTION:</b> Make sure that no parts of the tip funnel remain inside or on the catheter.</p>
18	<p>Press the THV entirely into the inflow cover tip.</p> <p><b>NOTE:</b> Avoid torquing or tilting the tip, only press the tip and bring it over the THV. <b>NOTE:</b> Make sure that the THV touches the catheter tip, otherwise the reference with the radiopaque ring will be wrong.</p>
19	<p>Completely close the catheter in the cold bowl until the edge of the cartridge and the tip touch. Remove the loading clip from the end of the delivery system. Move the switch on the handle from the position N (Neutral) to D (Drive). Repeatedly squeeze the lever of the handle until the handle touches the safety locker. Control the position of the THV. Move the switch on the handle from the position D (Drive) to N (Neutral) and move the handle forward until the flushing port of the catheter shaft is accessible (approximately 1 cm).</p> <p><b>NOTE:</b> Perform the test in cold sterile 0.9% saline solution in order to prevent warming and expansion of the valve. If expansion occurs, the centering hub can move and can prevent correct opening of the catheter shaft.</p>
20	<p>Use a 10 ml syringe, filled with heparinized sterile 0.9% saline solution and flush again the flushing port of the catheter shaft. Doing this, hold down the front part of the catheter (cartridge) in the bowl of the cold saline solution. Rinse until no further visible bubbles appear. Alternate between closing the catheter shaft and flushing until the edge of the cartridge and the catheter tip form a smooth transition.</p>
21	<p>Push the lever of the handle down and move the switch from the position N (Neutral) to the position P (Park). Remove the loading tool by moving it off the catheter.</p>
22	<p>Use a 10 ml syringe filled with heparinized sterile 0.9% saline solution and flush the integrated sheath at the appropriate flushing port, which is located on the front part of the function unit, until saline solution comes out at the end of the integrated sheath.</p>

Procedure	
23	<p>Use a 10 ml syringe, filled with heparinized sterile 0.9% saline solution and flush the guidewire lumen at the appropriate flushing port, which is located on the back part of the delivery system, until saline solution comes out at the tip of the catheter.</p> <p><b>NOTE:</b> Make sure, that the safety wings are outside of the traction tube.</p> <p><b>NOTE:</b> Make sure not to push the tip actuator during flushing!</p> <p><b>CAUTION:</b> Keep the THV hydrated and de-aired by flushing the port of the catheter shaft until ready for implantation.</p>
24	<p>Activate the hydrophilic coating of the cartridge and catheter shaft by wiping the surface with a swab soaked in a sterile 0.9% saline solution.</p>

### 7.9. THV DELIVERY

Before starting the procedure the implant site shall be measured to confirm the measurements previously performed during diagnostics.

**NOTE:** Use and monitor systemic anticoagulation during the implantation procedure based on hospital or physician protocol. If heparin is contraindicated, consider alternative anticoagulation therapy.

In the case of an implantation in a failing surgical bioprosthesis:

**CAUTION:** The type of the surgical bioprosthesis and the nominal diameter must be determined either by imaging or by reviewing patient records. Also determine the inner diameter and radiopaque marker of the degenerated surgical bioprosthesis using information from manufacturer information or literature.

**CAUTION:** Carefully check the patient's anatomy regarding low coronary artery location, especially in combination with shallow sinuses.

**CAUTION:** Ensure that there is no partially detached valve leaflet, which can occlude a coronary ostium in the aortic position.

**CAUTION:** The risk of coronary occlusion is higher for surgical bioprostheses with external leaflets.

**CAUTION:** Ensure that the guidewire is in the central lumen of the failing surgical bioprosthesis, because paravalvular placement or placement through a defect in the leaflet or the suture ring of the failing surgical bioprosthesis will lead to complications.

Procedure	
1	<p>The patient requires a contra and ipsilateral port with an introducer sheath in one side and a device compatible sheath at the other. Prepare the vascular access site according to standard technique.</p> <p><b>CAUTION:</b> Ensure the diameter of the access vessels is adequate.</p>
2	<p>Introduce a pigtail catheter under fluoroscopic guidance through the compatible sheath and position it behind the right or non-coronary leaflet.</p>
3	<p>Catheterize the left ventricle using an angiographic catheter in combination with a 0.035 inch (0.889 mm) guidewire with soft tip. Cross the aortic valve per preferred technique.</p>

Procedure	
4	<p>Perform an exchanging maneuver, remove the guidewire with the soft tip and advance the recommended 300 cm extra stiff guidewire through the catheter into the left ventricle. Pre-form the guidewire tip to an atraumatic J-shape prior insertion.</p> <p><b>CAUTION:</b> Make sure that if you pre-bend the guidewire it will not be damaged</p>
5	<p>Remove the catheter, leaving the guidewire in place in the left ventricle, and advance the valvuloplasty catheter with an appropriate diameter over the guidewire, cross the aortic valve, and position the balloon within the annular ring.</p> <p><b>NOTE:</b> The valvuloplasty catheter should be sized at least 1 to 2 mm smaller than the annular diameter. For the ALLEGRA 31 at least a 24 mm balloon is recommended.</p>
6	<p>Predilate the native aortic valve by balloon expansion under a short episode of rapid ventricular pacing. If a balloon valvuloplasty is performed within a failing surgical bioprosthesis, the diameter of the balloon should be less than the inner diameter of the failing surgical bioprosthesis. In addition, it is important to ensure that the balloon catheter is in the central lumen of the failed surgical bioprosthesis.</p> <p><b>CAUTION:</b> Balloon valvuloplasty should only be performed when the bioprosthesis is loaded into the delivery system and ready for implantation.</p> <p><b>CAUTION:</b> In case of failing surgical bioprostheses with structural damage, balloon valvuloplasty should be performed with extreme caution.</p> <p><b>NOTE:</b> The guidelines of the European Society of Cardiology state that balloon valvuloplasty should be contraindicated in left-sided degenerated surgical bioprostheses. Literature reports occurrences of fractures of surgical bioprosthetic support structures after balloon valvuloplasty.</p>
7	<p>Perform an aortic root angiography and find the best view of the valve to acquire a position where the three leaflets (NCL, RCL &amp; LCL) are aligned. When implanting in a failing aortic valve bioprosthesis, the radiopaque markers must also be used as guidance.</p>
8	<p>Introduce the delivery system, flushed, with the selected THV loaded on it and with the handle switch in the position P, over the guidewire under fluoroscopic guidance, throughout the introducer sheath to the aorta while maintaining guidewire position across the aortic valve.</p> <p><b>NOTE:</b> Before introducing the delivery system into the patient, visually check the correct loading of the THV.</p>
9	<p>Introduce the delivery system until the THV crosses the aortic valve.</p> <p><b>CAUTION:</b> If you feel any resistance, the procedure must be stopped and the cause of the resistance evaluated.</p>

## 7.10. THV DEPLOYMENT

Procedure	
1	Position the catheter so that the radiopaque ring of the inflow cover tip is on one level with the valve annulus marked by the pigtail catheter. In the case of an implantation in a failing bioprosthesis, position the catheter to align the THV with the position of the failed surgical bioprosthesis. The height of the implantation should be such that the gold markers are located at the distal end of the failing surgical bioprosthesis and the inflow portion of the THV is aligned with the inflow portion of the failing surgical bioprosthesis.
2	<p>Step 1: Opening</p> <p>With the switch in position D repeatedly squeeze the lever of the handle until the handle touches the safety locker.</p> <p>The catheter shaft will open and the THV will expand its middle portion and remain closed at its proximal and distal end (Permaflow).</p> <p><b>NOTE:</b> In case the catheter shaft has not fully opened, perform gentle push-pull movements until the catheter has fully opened.</p> <p><b>CAUTION:</b> Make sure that the catheter at the integrated sheath is stable, once positioned in correct annular position in order to avoid any change.</p>
3	Perform gentle push-pull movements of the catheter to precisely position and align the THV to the annulus plane. Position the radiopaque ring, located at the inflow cover tip, and align with the lowest (most proximal) side of the pigtail catheter loop which was positioned behind the native leaflets. In the case of an implantation in a failing bioprosthesis, perform gentle push-pull movements of the catheter to precisely align the THV to the desired implantation height within the failing surgical bioprosthesis.
4	<p>The six radiopaque markers of the stent at the level of the valve occluder will be in an expanded diameter and will be the reference to the position of the coronary ostium.</p> <p><b>NOTE:</b> The upper limit of the THV skirt is marked by the six radiopaque markers and must not be placed above the coronary arteries.</p>
5	The movement of the catheter to reposition the THV will also modify the angle of the THV in relation to the site of implantation. After the repositioning of the THV, gently move the catheter until the THV is realigned to the site of implantation.
6	<p>Step 2: Release</p> <p>After the correct alignment of the THV to the implantation site and confirmation that the coronary ostium is patent, securely hold the catheter at the integrated sheath to keep the position of the THV and move forward the tip actuator against the tip grip until these touch each other. The THV will be released and the proximal inflow end opens. The valve will immediately be functional.</p> <p><b>CAUTION:</b> Keep the delivery system very stable, in order not to dislocate the THV.</p>

Procedure	
7	In this phase the distal end of the THV with the T-Bars remains securely attached to the catheter. During the release of the tip actuator, the catheter tip will move inside the THV. <b>NOTE:</b> In case the THV must be retrieved, this can still be done after the release of the proximal inflow end. Follow the steps of safe bailout procedure.
8	Only when the tip comes inside the THV, the safety wings will retract and allow the release of the safety locker. <b>NOTE:</b> If for any reason the tip does not come inside the THV, the safety wings will not allow the release of the safety locker and the deployment should not be continued until this issue is solved. In this case repeat the actuation of the tip. <b>NOTE:</b> Make sure that the guidewire does not hinder the withdrawal of the catheter tip.
9	Gently pull the pigtail catheter to a position distal to the THV within the ascending aorta and perform an injection of contrast to confirm the position of the valve and the patency of the coronaries. <b>NOTE:</b> This must be performed under fluoroscopic control and making sure not to move the THV.
10	Unlock the safety locker after the clear assessment of the valve's correct position, by pressing the safety locker triggers simultaneously inward and moving the safety locker backward.
11	Step 3: Detachment Repeatedly squeeze the lever of the handle to withdraw the catheter shaft until the end of its course to detach the THV from the catheter totally. <b>CAUTION:</b> Check if the T-bars are free from the catheter. <b>NOTE:</b> Always pay attention to the stability of the THV.
12	Observe under fluoroscopy that the radiopaque ring, located at the cartridge end, has completely crossed the loading connector, which is also radiopaque.

### 7.11. WITHDRAWAL OF THE DELIVERY SYSTEM

Procedure	
1	Slowly withdraw the delivery system and observe under fluoroscopy if the delivery system tip is not trapped or dragging the THV, mainly when crossing the T-Bars.
2	Retract the guidewire, bringing its soft tip closer to the implanted valve to make this maneuver easy. This will take away the spring force placed by the guidewire against the catheter tip and therefore make the crossing of the valve safer.

Procedure	
3	<p>After the tip crosses the implanted valve, gently retract the delivery system into the descending aorta and close the catheter to be withdrawn. Set the handle switch to the neutral position N and hold the front cap of the delivery system with one hand. Slowly move the handle forward with the other hand. Advance the handle until the radiopaque marker band of the cartridge and the marker band of the inflow cover tip form a line but do not overlap the radiopaque tip.</p> <p><b>CAUTION:</b> Do not force the catheter against the tip to avoid overlapping. If this occurs the diameter of the tip will be increased by the overlapping catheter and make it difficult to cross the introducer sheath.</p> <p><b>NOTE:</b> The catheter closure is recommended to be done in a straight portion of the descending aorta.</p>
4	<p><b>NOTE:</b> The closure allows for withdrawal of the delivery system out of the introducer sheath without excessive bleeding.</p>

## 7.12. VERIFICATION OF PROSTHETIC VALVE POSITION AND MEASUREMENTS

Procedure	
1	<p>Perform a supra aortic angiogram to evaluate device position, function and coronary patency.</p> <p><b>CAUTION:</b> Proceed with caution if it is necessary to post-dilate the THV in the case of unacceptable paravalvular leakage. Post-dilatation could damage the device integrity or cause migration of the THV.</p> <p><b>NOTE:</b> If there is doubt about the coronary patency, an additional angiogram of each coronary branch can be selectively made using an angiographic catheter.</p>
2	<p>Measure and record the transvalvular pressure gradients.</p>
3	<p>Remove all catheters and sheaths.</p>
4	<p>Perform the closure of the arterial access according to preferred techniques and devices.</p>

## 7.13. SAFE VALVE RETRIEVAL PROCEDURE

In case the THV must be recaptured and retrieved from the patient, this can only be done after the release of the proximal inflow end and before further implantation steps are performed.

**CAUTION:** Re-sheath the THV only in bailout situations. Do not re-implant the THV after re-sheathing.

Procedure	
1	<p>Actuate the tip release to move the tip forward and fixate the position by using the loading clamp.</p>
2	<p>Move the switch from D to N and hold firmly the system.</p>
3	<p>To re-sheath move continuously the handle forward to recapture the THV into the cartridge until resistance occurs (approx. 1-2cm).</p> <p><b>CAUTION:</b> Apply only gentle force to avoid damage to the catheter. Keep the force applied until the end of the procedure.</p>

Procedure	
4	Carefully pull back the catheter through the aortic arch and stop if resistance occurs to verify its cause. Use the introducer sheath to complete the re-sheath of the THV.

## 8. HOW SUPPLIED

### 8.1. AVAILABLE SIZES

[mm]	ALLEGRA 23	ALLEGRA 27	ALLEGRA 31
Frame height	37.3	41.3	43.0
Inflow diameter	23.8	27.4	31.0
Outflow diameter	20.8	24.0	24.0
Aortic annulus diameter	19 - 22	22 - 25	25 - 28

### 8.2. PACKAGING

The THV is sterilized with and stored in glutaraldehyde solution. The THV is supplied sterile in a sealed product jar with a screw cap. The outside of the product jar is non-sterile and must not be placed in the sterile field. A warm and a cold temperature indicator are placed inside the carton box. **NOTE:** The product has been exposed to temperatures above 40°C if indicator is black. The product has been exposed to temperatures below 0°C if indicator bulb is purple. **CAUTION:** Do not use the THV if the warm temperature indicator is black or the cold temperature indicator's bulb is purple. **NOTE:** After implantation of the THV, the glutaraldehyde solution should be disposed of according to hospital procedure.

The delivery system and the loading system are sterilized with ethylene oxide gas (ETO). They are packed together in a double sterile bag. The delivery system and the loading system are sterile if the sterile bags are undamaged and unopened and the ETO indicator is green. The outer surface of the outer bag is non-sterile and must not be placed in the sterile field. **CAUTION:** Use the delivery system only if the ETO indicator on the inner sterile bag is green. **NOTE:** The used delivery system may be disposed of in the same manner that the hospital waste and bio-hazardous materials are handled. There are no special or unusual risks related to the disposal of the delivery system.

### 8.3. STORAGE

The THV must be stored between 10°C and 38°C. Avoid locations where extreme temperature fluctuation may occur. **CAUTION:** Do not freeze. Always store the THV in a cool, dry place. Any THV that has been frozen must not be used for human implantation.

The delivery system should be stored in a cool, dry place between 10°C and 38°C.

## 9. MRI SAFETY INFORMATION

Nonclinical testing indicates that the ALLEGRA Transcatheter Heart Valve is MR Conditional. It can be scanned under the following conditions:

- Static magnetic field of  $\leq 3$  Tesla
- Spatial gradient field of  $\leq 1500$  Gauss/cm

- Maximum whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning.
- Normal mode of MR system

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the THV.

#### 10. INFORMATION ON THE MATERIALS AND SUBSTANCES TO WHICH THE PATIENT CAN BE EXPOSED

Components and materials of the THV to which the patient will be exposed after implantation:

Component	Material	Quantitative composition
<b>Leaflets, Skirt</b>	Bovine pericardium (cross linked with glutaraldehyde)*	100%
<b>Suture material</b>	Polytetrafluoroethylene (PTFE)	100%
<b>Support stent structure</b>	Nitinol	100%
<b>Radiopaque markers</b>	Gold	100%

\* As the device contains animal tissue which is rendered non-viable, it is subject to Commission Regulation (EU) No 722/2012.

Components and materials of the Delivery System:

Component	Material
<b>Handle</b>	Polycarbonate (PC)
<b>Catheter</b>	Polyethylene (PE), Polycarbonate (PC), Polyetheretherketone (PEEK), Polyamide (PA), Polyimide (PI), Polytetrafluoroethylene (PTFE), Polyetherblockamide (PEBAX), Polyurethane (PU), Polyoxymethylene (POM)
<b>Hydrophilic coating</b>	Hydromer® Coating Solution S2005B
<b>Radiopaque markers</b>	Platinum/Iridium
<b>Others</b>	GoSilicone, Stainless steel

#### 11. PATIENT INFORMATION

An implant card and patient leaflet for the patient is included in each product package of the implantable device. The implant card contains blank fields to enter the patient name, date of implantation and healthcare institution filled by handwritten text. Information to identify the implanted device including device name, device model (device type), serial number and UDI as well as manufacturer address and website is provided by suitable implant card device stickers. Stickers are included within the device package of the implantable device. The implant card contains a dedicated field to place the implant card device sticker. The patient leaflet includes information regarding warnings, precautions or measures, expected lifetime of the device and necessary follow-up information. It also contains information for a healthcare professional with

regard to reciprocal interference with reasonably foreseeable external influences, medical examinations or environmental conditions. In addition, the overall qualitative and quantitative information on the materials and substances to which patients can be exposed is included in the patient leaflet. **NOTE:** The healthcare institution must ensure that the patient receives an implant card completely filled with all the predefined information as well as the information given by the patient leaflet. **NOTE:** Only a complete filled implant card allows the patient to have access to device related information.

## 12. TRAINING

NVT provides detailed loading and operator training. For training material and in-service support, contact a NVT representative. **CAUTION:** Implantation of the ALLEGRA Transcatheter Heart Valve may only be performed by physicians who have received the appropriate training.

**CAUTION:** Loading of the ALLEGRA Transcatheter Heart Valve may only be performed by physicians or hospital personnel who have received the appropriate training.

## 13. REPORTING TO MANUFACTURER AND COMPETENT AUTHORITY

Any health care professional (e.g. customer or user of the device) and/or a patient treated with the ALLEGRA THV who has identified any serious incident that has occurred in relation to the provided device should inform the manufacturer and the local competent authority immediately. The manufacturer has to be notified by telephone, FAX, or in written correspondence. When reporting a serious incident provide device name and model, lot number or serial number and the nature of incident.

## 14. SUMMARY OF SAFETY AND CLINICAL PERFORMANCE AND EXPECTED CLINICAL BENEFIT

For implantable devices, the manufacturer sets up a regular report of safety and clinical performance (SSCP). A summary of this report for healthcare professionals and patients will be made publicly available in the future via Eudamed (European Databank on Medical Devices).

Until the Eudamed is fully functional the SSCP can be requested at the website of the manufacturer. The clinical benefit of the ALLEGRA Transcatheter Heart Valve has been proven by multiple pre-clinical trials to show the safety and efficacy according to the intended use of the device. Clinical results were assessed by medical experts and notified body, which have led to the CE-mark approval followed by the commercial use of the device. On-going studies are conducted to show the clinical performance of the device. The results are summarized and updated in the regular report of safety and clinical performance (SSCP).

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Rev.00070301\_082023\_13\_ML

### Disclaimer

ALLEGRA is a trademark or registered trademark of NVT AG. ALLEGRA™ is CE Mark approved. CAUTION: The law restricts these devices to sale by or on the order of a physician and these products are intended for the use by or under the direction of a physician. Not available in the United States and any other country where applicable health authority product registration has not been obtained. Information contained herein only for presentation outside the US and France. ALLEGRA™ is a product of NVT AG. Blue Sail Medical Co., Ltd is the parent company of NVT AG and Biosensors Intl. and its affiliates are collaborating for the commercialization of the ALLEGRA™ device. © 2021 NVT AG. All rights reserved.

### NVT GmbH

Lotzenäcker 17  
72379 Hechingen (Germany)  
Tel. +49 (0) 7471 98979-0  
info@nvt-med.com

[www.nvt-med.com](http://www.nvt-med.com)