en: Patient Leaflet

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General:

Note: Please read all information carefully provided by this patient leaflet. The manufacturer provides with each implant an implant card and a patient leaflet. A completely filled implant card provides important information about the implant, which has been used for the treatment. The filled implant card must be stored carefully by the patient and allows safe traceability of the implant for the patient. The patient leaflet includes information regarding warnings, precautions or measures, expected lifetime of the device and necessary follow-up information. It also contains information for a healthcare professional with regard to reciprocal interference with reasonably foreseeable external influences, medical examinations or environmental conditions. In addition, the overall qualitative and quantitative information on the materials and substances to which patients can be exposed is included.

Information for the patient:

Implant Card: Make sure that the implant card is completely filled with all required data when provided by the health care institution. Keep the implant card and the patient leaflet carefully and show to the healthcare institution when requested during follow-up examination.

Note: Keep the implant card and patient leaflet carefully in your personal medical documentation.

Expected life time: At some time the ALLEGRA transcatheter heart valve may need to be replaced. Careful and continuous medical follow-up is advised so that implant related complications can be diagnosed and properly managed. The ALLEGRA transcatheter heart valve was tested to mimic 5 years of use and even tested for 10 years (metal stent frame) according to current applicable standards.

Follow-up information: Regular follow-up examination are important. You will probably be asked to return to your responsible healthcare institution to have your transcatheter heart valve checked after your procedure (often at 30 days and 1 year). Follow-up examination may vary.

Note: Contact your responsible healthcare institutions whenever you have questions, concerns or changes in your health status.

Precaution: It is important to share information about your transcatheter heart valve replacement before any medical, dental, or MRI (magnetic resonance imaging) procedures or examinations. If you need an MRI provide information that you have received an ALLEGRA transcatheter heart valve and the MRI information given in this patient leaflet.

Precaution: Make sure that you have received and discussed your aftercare plan with the physician. The health-care institution may provide you with specific instructions to help you with your recovery. This may include a special diet, when to return to exercise, and any medicine you may need to take.

Note: It is important to carefully follow the directions given by the physician, especially if you need to take any blood thinning medication. Don't change the medication yourself without consulting your doctor.

Warning: There are no specific warnings related to the ALLEGRA transcatheter heart valve itself. However, consider all additional warnings which may be provided by your healthcare institution.

Information for the healthcare professional:

Implant card: Fill in the name of the patient. Add implantation date and the address of the responsible physician and/or your healthcare institution. After treatment, place the implant card device sticker on the dedicated place on the backside of the implant card. The implant card device sticker can be found in the product package of each ALLEGRA transcatheter heart valve. After completion, handover the implant card and the patient leaflet to the patient.

Note: Make sure that all information required by the implant card is correctly and readable filled and complete.

Warning: It is recommended to the treating institution that all transcatheter heart valve recipients be prophylactically treated for endocarditis according guidelines to minimize the risk of prosthetic valve infection.

Reasonable foreseeable external influences or environmental conditions: The ALLEGRA transcatheter heart valve is MR conditional. MRI examination can be performed according to the given specification in this document.

Medical examination: Consider the ESC guideline regarding management of infective endocarditis. (2015 ESC Guidelines for the management of infective endocarditis, European Heart Journal, 2015)

Reciprocal interference: There are no known reciprocal interference in relation to the ALLEGRA transcatheter heart valve.

Materials and quantitative composition:

Materials and quantitative composition to which the patient will be exposed after device implantation:

Component	Material	Quantitative composition
Leaflets, Skirt	Bovine pericardium	100%
Suture material	Polytetrafluoroethylene (PTFE)	100%
Support stent structure	Nitinol	100%
Radiopaque markers	Gold	100%

MRI Safety Information:

Non-clinical testing indicates that the ALLEGRA transcatheter heart valve is MR Conditional. A patient with these implant can be scanned safely under the following conditions: Static magnetic field of ≤ 3 Tesla, Spatial gradient field of ≤ 1500 Gauss/cm, Maximum whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning. Normal mode of MR system was chosen. MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the implant.

Manufacturer Information:

This patient leaflet was revised the last time in March 2020. Regular update information will be provided at company homepage (www.nvt-med.com.) under the chapter patient information.

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Copyright: NVT GmbH. All rights reserved.
NVT GmbH